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Data are presented as the means ± SD. GFR, glomerular filtration rate; NAG, N-acetyl-β-D-glucosaminidase; β2-MG, β2 microglobulins; L-FABP, liver-type fatty acid-binding protein; CI-AKI, contrast-induced acute kidney injury. Non-CI-AKI group vs. CI-AKI group: \*P<0.05; \*\*P<0.001.

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- Komatsu S, Hirayama A, et al. Detection of Coronary Plaque by Computed Tomography with the Novel Plaque Analysis System "Plaque Map" and Comparison with Intravascular Ultrasound and Angioscopy. *Circ J* 2005, 69 (1), 72-77.
- 小松誠:緊急冠動脈CTの臨床的意義 児玉和久監修、平山篤志編集、Stephan Achenbach協力編集「循環器診療・インターベンションのためのMDCT」p83-87 中山書店 2007.
- Cigarroa RG, Lange RA, Williams RH, Hillis LD. Dosing of contrast material to prevent contrast nephropathy in patients with renal disease. *Am J Med.* 1989 Jun;86(Pt 1):649-52.
- Komatsu S, Kamata T, Imai A, Ohara T, Takewa M, Miyaji K, Yoshida J, Kodama K. Controlling Intracoronary CT Number for Coronary CT Angiography. *J Cardiol* 2013 Feb;61(2):155-61
- Komatsu S, Kamata T, Imai A, Ohara T, Takewa M, Ohe R, Yoshida J, Kodama K. Coronary computed tomography angiography using ultra-low-dose contrast media: radiation dose and image quality. *Int J Cardiovasc Imaging* 2013, 29, 1335-1340.
- Veerkamp JH, Peeters RA, Maatman RG. Structural and functional features of different types of cytoplasmic fatty acid-binding proteins. *Biochim Biophys Acta.* 1991 Jan 4;1081(1):1-24.
- Yokoyama T, Kamijo-Ikemori A, Sugaya T, Hoshino S, Yasuda T, Kimura K. Urinary excretion of liver type fatty acid binding protein accurately reflects the degree of tubulointerstitial damage. *Am J Pathol.* 2009 Jun;174(6):2096-106.
- Matsui K, Kamijo-Ikemori A, Sugaya T, Yasuda T, Kimura K. Usefulness of urinary biomarkers in early detection of acute kidney injury after cardiac surgery in adults. *Circ J.* 2012;76(1):213-20.
- Manabe K, Kamihata H, Motohiro M, Senoo T, Yoshida S, Iwasaka T. Urinary liver-type fatty acid-binding protein level as a predictive biomarker of contrast-induced acute kidney injury. *Eur J Clin Invest.* 2012 May;42(5):557-63.
- Komatsu S, Takahashi S, Kodama K: Urinary Liver Fatty-acid Binding Protein within 24 Hours from Contrast Media Exposure May Predict Acute Kidney Injury. American Heart Association Scientific Sessions, Nov 19, 2014 (Chicago).